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Product features:

A) Description:

Drills and instruments for dental implants.

B) Materials:

Items are made of 17 4PH (AISI 630) hardened stainless steel, electro-sanitized with specific acid bath and cleaned with ultrasound, or of Widia K 30 with physical/mechanical specifications specific for use according to current regulations.

C) Instructions for use:

Burs and extensions must be used with the help of a handpiece (micro motor). Recommended working speed is from several rpm to 800 rpm (high speed use over 1000 rpm is discouraged), except for the bone tap and screwdrivers where the speed must be set to the minimum available of your handpiece or in any case not higher than 10 rpm.

D) Packaging features:

Items are packaged in bags suited to preserve functional device characteristics and cleanliness.

Packaging symbol legend:

Manufacturer

Do not use if package is damaged

Consult instructions for use

Catalogue number

Non-sterile

Batch code

CE Mark

0476

E) Warnings:

- Extension use is recommended with SAEG burs or with burs that meet UNI EN ISO 1797-1 regulation tolerance.
- All instruments must be cleaned (ultrasound equipment is recommended) and sterilized before use as per regulation UNI EN 556-1.
- Cleaning and sterilization must only be performed with specific materials, specifically, detergents containing the following substances should be avoided:
 - Oxalic acid
 - High concentrate chlorine.
- Insert the instruments in the disinfectant solution immediately after surgery to prevent encrustation of blood, secretions, etc.
- Do not store wet or damp instruments.
- Do not sterilize, clean or disinfect instruments made with different metals in the same cycle. Steam sterilization in autoclave is specifically recommended for cutting instruments to prevent dulling.
- Al In order to guarantee optimal bur performance, it should not be used more than ten times.
- Diameter notches are included where permitted by piece dimensions and physiognomy.

F) Product identification:

SAEG products are identified by the product code and production lot printed on the packaging, while on the device they are marked (where the dimensions allow it) diameter and length of the drill. One type of drills also has a color band for faster identification of the diameter.

G) Specific liability:

SAEG oral implant products meet the essential requisites set by directive 93/42/CEE and directive 2007/47/CE. Any product use other than that specified is considered "improper" and releases the manufacturer from any liability.

H) Possible causes for product malfunction:

- The bur resists entrance:

Check bur sharpness, it may be dull due to repeated use and sterilization.

- The bur overheats upon entrance:

Make sure internal irrigation with physiological solution is adequate to avoid contacted cell necrosis.

- The bur is broken:

Remove all fragments from the site and adequately clean the incision before continuing..

- The bur does not lock onto the extension:

Check bur surface and extension surface tolerance which may be worn due to overuse.

- The interchangeable stop does not remain inserted in the special seat on the drill:

The centesimal tolerance may have been compromised by incorrect maintenance / cleaning or by repeated use. Replace the device with a new one.

WARNING	Saeg's cutting instruments need special attention during the cleaning to do not damage the sharpening.
Restriction	The end of the tool life has brought about normal wear and damages for the uses.
INSTRUCTION	
Use way	Devices should be cleaned and sterilized before every use.
Box and transport	The used box for the transport guarantee the correct cleanliness for the sterilization but they are not suitable for it.. If the cleanliness and sterilization process are made outside own company is better to use the plastic box for the transport
Cleanliness preparation	All the devices with moving parts must be disassembled, see specific use instruction. It is better and recommended to make the preliminary cleanliness and the decontamination by ultrasound apparatus, 20 minutes 70°C temperature using suitable detergent for stainless steel medical devices.
Cleaning by machine	Make the cleaning cycle with instrument washing machine respecting the instruction of the machine producer, checking very well the instrument position.
Manual cleaning	Remove by soft brush the bone/blood residual checking very well the hole for internal cooling.
Disinfection	Put the cleaned instruments into BIOSAN UNO solution for 3 minutes (see the instruction of use).
Drying	Not necessary (eventually by compressed air).
Maintenance	Verify the drill cut if it is not suitable change the one.
Inspection and functional tests	Check that the cutting part of the device does not show any damage and signs of wear. The cutting edges must be uniform and free from chipping. Check that any internal irrigation hole (with the respective outlet hole) is not blocked by any residues. Verify that the assembly systems are still suitable.
Packaging	Singularly: check the packaging material be enough strong and the drill cut do not damage it.. More pieces: to use special trays for instruments or sterilization ones for general use. To wrap the tray using the suitable way.
Sterilization	By autoclave for 15 minutes minimum (134°C).
Storage	Clean and dry place.
Added information	For the sterilization by only cycle, verify to have not exceed the maximum charge of the sterilizer.
Manufactured contact	Very important is to know the phone number of the manufacture to solve possible problems.

The instructions provided above have been validated by the medical device manufacturer as ABLE to prepare a medical device for reuse. The process manager has the responsibility to ensure that repeated processes are actually performed using the equipment, materials and personnel in the repeat process facility to achieve the desired result. This generally requires systematic validation and monitoring of the process. Similarly, all deviations by the process manager from the instructions provided should be adequately assessed to judge their effectiveness and potential undesirable consequences.

Incorrect or insufficient maintenance can damage instruments in a short time.

There are specific products on the market for cleaning instruments, strictly observe the instructions for use.

**WARNING: SAEG DEVICE USE IS RESERVED TO DENTAL SURGEONS.
THE CHOICE OF THE DEVICE MUST BE CONSEQUENTIAL TO A ACCURATE HISTORY OF THE PATIENT.**

For information, clarifications or suggestions, please contact SAEG.

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