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Manufacturer Identification:

The Manufacturer and Person Responsible for placing on the EU territory the Medical Devices referred to as DRILLS AND DENTAL SURGICAL INSTRUMENTS covered by this Instruction for Use is:

SAEG S.r.l.
C.so G. Di Vittorio 5127029 Vigevano (PV) ITALY
Tel. +39 0381 76350 - Fax +39 0381 72667
www.saeg.it info@saeg.it

Device Identification:

Lance Drill, Pilot Drill, Twist tri spade, DLC-Coated Drill, TIN-Coated Drill, Drill Extensions, Drill Stop, Threaded Drill Stop, Hand tissue punch, Crestal Pin, Countersink, Core Drill, Screw Driver and Tap Handpiece.

Intended use:

The DRILLS AND DENTAL SURGICAL INSTRUMENTS are intended for the preparation of the surgical socket in which the implant is to be inserted, allow the extraction of the implant from the package, its transport to the patient's mouth and the insertion of the implant itself, it also allows the correct positioning and insertion of the prosthetic components on the implant

SAEG products are identified with the product code and production batch on the packaging, while the diameter and drill length are marked on the device (where size permits). If requested by the customer, SAEG places a colour band on the drills for faster identification of the diameter.

Materials:

Items are made of 17 4PH (AISI 630) hardened stainless steel, electro-cleaned with a specific acid bath and subsequent ultrasonic cleaning, with physical-mechanical characteristics suitable for use, in accordance with current regulations.

Method of use:

Drills, taps and extensions must be used with the aid of a handpiece (micromotor) the recommended working speed is from a few revolutions to 800 rpm (use at speeds above 1000 rpm is not recommended) except for handpiece taps and screwdrivers where the speed must be set to the minimum available of your handpiece or in any case not more than 10 rpm.

Maximum number of uses of the devices N. 10.

Packaging features:

The items are packaged in bags designed to maintain the appropriate degree of cleanliness for proper sterilisation

Warning:

We recommend using the extension with SAEG drills or with drills that comply with the tolerances of UNI EN ISO 1797-1.

All instruments must be cleaned (ultrasonic equipment is recommended) and sterilised before use, following the instructions of the steam steriliser manufacturer.

Cleaning and sterilisation must only be carried out with specific materials. Do not use cleaning products containing the following:

- Oxalic Acid
- Chlorine in high concentrations.
- Place the instruments in the disinfectant solution immediately after surgery and leave them for a few hours to prevent any encrustation of blood, secretions etc.
- Do not store wet or damp instruments.
- Do not sterilise, clean or disinfect instruments made with different metals in the same cycle.

For sharp instruments, steam sterilisation in an autoclave is particularly recommended so as not to deteriorate the sharpness of the cutting parts.

In order to ensure optimal functionality of the drills, we recommend no more than ten uses.

Diameter marking is carried out where the dimensions and physiognomy of the workpiece allow it.

Users:

The use and handling of the medical device is reserved for medical and dental personnel with the necessary qualifications and professional training.

THE CHOICE OF THE DEVICE MUST FOLLOW A THOROUGH REVIEW OF THE PATIENT'S MEDICAL HISTORY

Specific liability:

SAEG products for oral implantology are compliant with the GENERAL SAFETY AND PERFORMANCE REQUIREMENTS of Regulation (EU) 2017/745-MDR. Any use of the products other than that specified is to be considered as "improper use" and relieves the manufacturer of any liability.

POSSIBLE CAUSES OF PRODUCT MALFUNCTION:

The drill forces its way into the site

Check the sharpness of the drill, it may have lost its cutting characteristics due to repeated use and sterilisation. Replace the device with a new one.

The drill entering the site heats up excessively

Check that internal or external irrigation with physiological substance is adequate to avoid necrosis of the insulted cells.

The drill has broken

Remove all fragments from the site and clean the wound adequately before continuing the procedure with a new drill.

The drill does not remain locked in the extension

Check the tolerance of the drill table and that of the extension, which may have become worn due to too many uses.

The interchangeable stop does not remain inserted in the seat on the drill

The centesimal tolerance may have been compromised by incorrect maintenance/cleaning or by repeated use. Replace the device with a new one.

INSTRUCTIONS FOR CORRECT USE AND MAINTENANCE

WARNINGS	Saeg cutting tools require special care when cleaning to avoid damaging the sharpness. Internally irrigated drills require special care when cleaning.
Limitation to repeated processes	The end of the service life is generally determined by wear and damage due to use; however, incorrect maintenance may alter the functionality of the device.
Site of use	The devices supplied must be cleaned and sterilised before each use.
Containers and transport	Containers used for transport provide the appropriate degree of cleanliness for proper sterilisation but are not suitable for sterilisation. If cleaning and sterilisation processes are carried out outside the facility, it is recommended that these containers be reused so as not to damage the devices during transport.

Preparation for cleaning	All devices with moving parts must be disassembled. Pre-cleaning and decontamination with ultrasonic equipment is recommended for a minimum of 20 minutes at a temperature not exceeding 70°, using a detergent suitable for cleaning stainless steel medical equipment.			
Cleaning: Automated	Carry out the cleaning cycle with a thermo-disinfecting washer according to the manufacturer's instructions and make sure that the instruments are loaded without shadowing (no overlapping).			
Cleaning: Manual	Remove bone/blood residues from the place of use with a soft brush. Check that the irrigation holes are not blocked and free them with a small needle if necessary.			
Disinfection	A disinfectant solution may be used for no longer than 3 minutes (as specified by the manufacturer), or Henry Schein generic detergent or equivalent product for cleaning dental surgical instruments			
Drying	Not required (possibly with a compressed air gun).			
Maintenance	Check that the cutting edge of the drill is suitable for further use and if not, discard the device.			
Inspection and function tests	Check the cutting edge of the device for damage and signs of wear. The cutting edges should be even and free of chips. Check that the internal irrigation hole (with its respective outlet holes), if any, is not blocked by residue. Check that the assembly systems are still suitable.			
Packaging	<p>Individually: a reference packaging material may be used. Make sure that the packaging is strong enough so that it will not be damaged by the sharp edge of the device.</p> <p>In groups: instruments can be loaded onto trays intended for instruments, or onto sterilisation trays for general use, without storing instruments of different materials. Make sure that sharp parts are protected. Wrap the trays using the appropriate method.</p>			
Sterilisation	Pressure	Time	Temperature	Drying
Cycle	2 Atm	3-7 minutes	134°C	20 minutes
Storage	In a clean and dry place.			
Additional information	Additional information For sterilisation in a single cycle, ensure that the maximum load of the steriliser is not exceeded.			
Manufacturer Contact	See brochure for telephone number and address of local representative.			

Disposal Procedures:

Surgical instruments, if removed from the oral cavity due to biological or mechanical failure, must be assimilated for disposal as biological waste, according to local regulations.

Reporting Accidents:

In accordance with SECTION 2 VIGILANCE, Article 87 Serious Accident Reporting and Safety Corrective Actions of MDR 745/2017/EU, Legislative Decree of 15 November 2005 and MEDDEV Guideline 2.12-1, it is recalled that public or private healthcare professionals are required to report accidents or near misses with Medical Devices to the manufacturer (SAEG S.r.l.) or to the National Competent Authorities as a matter of urgency.

An accident is defined as any malfunction or deterioration in the characteristics and/or performance of a device, as well as any deficiency in the labelling or operating instructions which, directly or indirectly, can cause or has caused death or serious deterioration of the state of health of the patient or a user or other person.

Further Information:

This device does not contain software.

This device must be sterilised before use (the device is supplied in NON-STERILE condition).

This device is intended to be reused.

This device is not intended for use by inexperienced users.

The device is not included in the list of product groups that are not intended for medical use.

Defective Product Liability And Warranty Terms:

Optimal patient care and attention to the patient's needs are necessary conditions for implant success, and it is, therefore, necessary to carefully select the patient, inform him or her of the inherent risks and duties associated with treatment, and encourage him or her to cooperate with the dentist for the successful outcome of the treatment. It is, therefore, necessary for the patient to maintain good hygiene, confirmed during check-ups and follow-up appointments; this must always be ensured and documented, just as the doctor's pre- and post-operative instructions and prescriptions must be observed and documented.


The instructions provided by SAEG S.r.l. are available at the time of treatment and accepted in dental practice; they must be observed and applied in all phases of treatment: from the patient's anamnesis to post-operative check-ups.

The warranty covers only ascertained production defects, after sending the part identified by part number and batch, within the warranty period. Warranty terms are available at www.saeg.com.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for re-use. It is the responsibility of the process manager to ensure that repeat processes are actually performed using the equipment, materials and personnel in the repeat process structure to achieve the desired result. This generally requires validation and systematic monitoring of the process. Similarly, all deviations by the process manager from the instructions provided should be properly evaluated to judge their effectiveness and potentially undesirable consequences.

Incorrect or insufficient maintenance can potentially damage instruments in a short period of time.

There are specific products on the market for cleaning these instruments. Make sure you follow the operating instructions correctly.

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Manufacturer's name		CE marking with notified body identification
	Date of manufacture		Production batch
	Consult the operating instructions		Attention! See the operating instructions
	Non-sterile device		Do not use if the packaging is damaged
	Device identification code		The product is a medical device
	Unique Device Identification - Alphanumeric code that uniquely identifies a medical device		

Date And Validity Of Operating Instructions:

These operating instructions are valid and effective from February 2023.

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