

Founded in 1973 as manufacturer of special tools, over the years Saeg has specialized in the production of medical devices for oral implantology.

Thanks to the experience in that field, Saeg offers to customers a high quality product.

The entire production is carried out by highly qualified operators using last generation CNC electronic machinery. All products made of stainless steel are subjected to a treatment called passivation which ensures a very high degree of resistance to oxidation.

Our production systems are certified ISO 9001 and ISO 13485. All the devices are manufactured following the latest revisions of the current regulation.



Last update:



October 2016



100%
MADE IN ITALY

SAEG STRENGTHS

- ▶ High quality of raw material
- ▶ Products made entirely by Saeg (except coatings)
- ▶ Realization of any drill and / or equipment from customer's drawing with possibility of CE marking
- ▶ Possibility of TIN and DLC coating
- ▶ Passivation of all steel devices
- ▶ Ability to laser engraving any kind of writing or customer's logo (where size permits)
- ▶ Possibility of color coding
- ▶ Ability to realize drills with internal irrigation
- ▶ Entirely made in Italy
- ▶ "Guided Implantology System" patented by Saeg



LANCE DRILL

Lance drill made by Hard Metal apt to open the cortical bone.

Article	Total Length	Useful Length
FIL-235	L 22.5	L 10
FIL-235C	L 26	L 13.5
FIL-235E	L 30	L 17.5
FIL-235L	L 37.5	L 25



PILOT DRILL

Pilot drill made by stainless steel apt to begin the implant site. Available with or without depth laser marks.

Article	Diameter	Total Length	Useful Length	Depth Marks
FIP-120	Ø1.2	30.0	15.5	On demand
FIP-150	Ø1.5	31.5	17	On demand
FIP-180	Ø1.8	34.5	22	On demand
FIP-200C	Ø2.0	28.5	16	On demand
FIP-200L	Ø2.0	34.5	22	On demand
FIP-225C	Ø2.25	28.5	16	On demand
FIP-225L	Ø2.25	34.5	22	On demand

Available with TIN or DLC coating.



TWIST TRI-SPADE

Twist tri-spade made by stainless steel apt to progressively enlarge the implant site until reaching the implant characteristics. Available with or without depth laser marks.

Article	Diameter	Total Length	Useful Length	Hole	Depth Marks
FIF-215D	Ø2.15	32	12.5	No	On demand
FIF-215M	Ø2.15	38	18.5	No	On demand
FIF-250D	Ø2.50	32	12.5	No	On demand
FIF-250M	Ø2.50	38	18.5	No	On demand
FIF-280D	Ø2.80	32	12.5	No	On demand
FIF-280M	Ø2.80	38	18.5	No	On demand
FIF-300D	Ø3.00	32	12.5	No	On demand
FIF-300M	Ø3.00	38	18.5	No	On demand
FIF-325D	Ø3.25	32	12.5	No	On demand
FIF-325M	Ø3.25	38	18.5	No	On demand
FIF-350D	Ø3.50	32	12.5	No	On demand
FIF-350M	Ø3.50	38	18.5	No	On demand
FIF-380D	Ø3.80	32	12.5	No	On demand
FIF-380M	Ø3.80	38 </td <td>18.5</td> <td>No</td> <td>On demand</td>	18.5	No	On demand
FIF-400D	Ø4.00	32	12.5	No	On demand
FIF-400M	Ø4.00	38	18.5	No	On demand
FIF-415D	Ø4.15	32	12.5	No	On demand
FIF-415M	Ø4.15	38	18.5	No	On demand
FIF-445D	Ø4.45	32	12.5	No	On demand
FIF-445M	Ø4.45	38	18.5	No	On demand
FIF-510D	Ø5.10	32	12.5	No	On demand
FIF-510M	Ø5.10	38	18.5	No	On demand
FIF-560D	Ø5.60	32	12.5	No	On demand
FIF-560M	Ø5.60	38	18.5	No	On demand



INTERNAL COOLING TWIST TRI-SPADE

Twist tri-spade made by stainless steel apt to progressively enlarge the implant site until reaching the implant characteristics. Made with internal irrigation holes. Available with or without depth laser marks.

Article	Diameter	Total Length	Useful Length	Hole	Depth Marks
FIF-215C	Ø2.15	32	12.5	Yes	On demand
FIF-215L	Ø2.15	38	18.5	Yes	On demand
FIF-250C	Ø2.50	32	12.5	Yes	On demand
FIF-250L	Ø2.50	38	18.5	Yes	On demand
FIF-280C	Ø2.80	32	12.5	Yes	On demand
FIF-280L	Ø2.80	38	18.5	Yes	On demand
FIF-300C	Ø3.00	32	12.5	Yes	On demand
FIF-300L	Ø3.00	38	18.5	Yes	On demand
FIF-325C	Ø3.25	32	12.5	Yes	On demand
FIF-325L	Ø3.25	38	18.5	Yes	On demand
FIF-350C	Ø3.50	32	12.5	Yes	On demand
FIF-350L	Ø3.50	38	18.5	Yes	On demand
FIF-380C	Ø3.80	32	12.5	Yes	On demand
FIF-380L	Ø3.80	38	18.5	Yes	On demand
FIF-400C	Ø4.00	32	12.5	Yes	On demand
FIF-400L	Ø4.00	38	18.5	Yes	On demand
FIF-415C	Ø4.15	32	12.5	Yes	On demand
FIF-415L	Ø4.15	38	18.5	Yes	On demand
FIF-445C	Ø4.45	32	12.5	Yes	On demand
FIF-445L	Ø4.45	38	18.5	Yes	On demand
FIF-510C	Ø5.10	32	12.5	Yes	On demand
FIF-510L	Ø5.10	38	18.5	Yes	On demand
FIF-560C	Ø5.60	32	12.5	Yes	On demand
FIF-560L	Ø5.60	38	18.5	Yes	On demand



CONTRA-ANGLE TISSUE PUNCH

Tissue punch made by stainless steel apt to create a circular incision of the gingival mucosa.

Article	Diameter	Total Length
CMU-300	Ø3.0	28
CMU-350	Ø3.5	28
CMU-410	Ø4.1	28
CMU-450	Ø4.5	28
CMU-510	Ø5.1	28

Available with TIN or DLC coating.



TREPHINE

Trephine made by stainless steel apt to remove implant and to take bone samples.

Article	Internal Diameter	External Diameter	Total Length	Depth Marks
FCR-300	Ø3.0	Ø3.7	33	5-10-15
FCR-400	Ø4.0	Ø4.7	33	5-10-15
FCR-500	Ø5.0	Ø5.7	33	5-10-15
FCR-600	Ø6.0	Ø6.7	33	5-10-15

Available with TIN or DLC coating.



CYLINDRICAL DRILL WITH DEPTH MARKS

Drill made by stainless steel with TIN or DLC coated 16mm useful length, with color-coded for diameter identification and laser depth marks on demand.

Article	Diameter	Color
FAB-340/01	Ø2.0	White
FAB-340/02	Ø2.5	Yellow
FAB-340/03	Ø2.8	Red
FAB-340/04	Ø3.2	Blue
FAB-340/05	Ø3.65	Green
FAB-340/06	Ø4.3	Black
FAB-340/07	Ø5.2	Violet



CYLINDRICAL DRILL WITH CALIBRATED LENGTH

Drill made by stainless steel available in 7 diameters for 5 lengths, with calibrated length.

Article	Diameter	Useful Length	Color
FAB-644/01	Ø2.0	6.0	White
FAB-644/02	Ø2.0	8.0	White
FAB-644/03	Ø2.0	10.0	White
FAB-644/04	Ø2.0	11.5	White
FAB-644/05	Ø2.0	13.0	White
FAB-644/06	Ø2.5	6.0	Yellow
FAB-644/07	Ø2.5	8.0	Yellow
FAB-644/08	Ø2.5	10.0	Yellow
FAB-644/09	Ø2.5	11.5	Yellow
FAB-644/10	Ø2.5	13.0	Yellow
FAB-644/11	Ø2.8	6.0	Red
FAB-644/12	Ø2.8	8.0	Red
FAB-644/13	Ø2.8	10.0	Red
FAB-644/14	Ø2.8	11.5	Red
FAB-644/15	Ø2.8	13.0	Red
FAB-644/16	Ø3.2	6.0	Blue
FAB-644/17	Ø3.2	8.0	Blue
FAB-644/18	Ø3.2	10.0	Blue
FAB-644/19	Ø3.2	11.5	Blue
FAB-644/20	Ø3.2	13.0	Blue
FAB-644/21	Ø3.65	6.0	Green
FAB-644/22	Ø3.65	8.0	Green
FAB-644/23	Ø3.65	10.0	Green
FAB-644/24	Ø3.65	11.5	Green
FAB-644/25	Ø3.65	13.0	Green
FAB-644/26	Ø4.3	6.0	Black
FAB-644/27	Ø4.3	8.0	Black
FAB-644/28	Ø4.3	10.0	Black
FAB-644/29	Ø4.3	11.5	Black
FAB-644/30	Ø4.3	13.0	Black
FAB-644/31	Ø5.2	6.0	White
FAB-644/32	Ø5.2	8.0	White
FAB-644/33	Ø5.2	10.0	White
FAB-644/34	Ø5.2	11.5	White
FAB-644/35	Ø5.2	13.0	White

All length are overestimated 0.3mm.



DRILL WITH INTERCHANGEABLE STOP

Drill made by stainless steel designed to put on interchangeable stop.

Article	Diameter
FSV-1/200	Ø2.0
FSV-1/250	Ø2.5
FSV-1/280	Ø2.8
FSV-1/300	Ø3.0
FSV-2/320	Ø3.2
FSV-2/350	Ø3.5
FSV-2/365	Ø3.65
FSV-2/400	Ø4.0
FSV-2/430	Ø4.3
FSV-3/450	Ø4.5
FSV-3/480	Ø4.8
FSV-3/520	Ø5.2
FSV-3/560	Ø5.6



Article	Length
SPF-1/08	8.0
SPF-1/09	9.0
SPF-1/10	10
SPF-1/11.5	11.5
SPF-1/12	12
SPF-1/13	13
SPF-1/14	14
SPF-2/08	8.0
SPF-2/09	9.0
SPF-2/10	10
SPF-2/11.5	11.5
SPF-2/12	12
SPF-2/13	13
SPF-2/14	14
SPF-3/08	8.0
SPF-3/09	9.0
SPF-3/10	10
SPF-3/11.5	11.5
SPF-3/12	12
SPF-3/13	13
SPF-3/14	14



The useful length of the drill with assembled stop is overestimated 0.5mm.

DRILL WITH THREADED INTERCHANGEABLE STOP

Drill made by stainless steel designed to put on threaded interchangeable stop.

Article	Diameter	M
FFS-1/200	Ø2.0	M3.5
FFS-1/250	Ø2.5	M3.5
FFS-1/280	Ø2.8	M3.5
FFS-2/300	Ø3.0	M4.5
FFS-2/320	Ø3.2	M4.5
FFS-2/350	Ø3.5	M4.5
FFS-2/365	Ø3.65	M4.5
FFS-3/400	Ø4.0	M6
FFS-3/430	Ø4.3	M6
FFS-3/480	Ø4.8	M6
FFS-3/520	Ø5.2	M6



Article	Length	M
STF-1/06	6.0	M3.5
STF-1/08	8.0	M3.5
STF-1/10	10	M3.5
STF-1/11.5	11.5	M3.5
STF-1/13	13	M3.5
STF-1/16	16	M3.5
STF-2/06	6.0	M4.5
STF-2/08	8.0	M4.5
STF-2/10	10	M4.5
STF-2/11.5	11.5	M4.5
STF-2/13	13	M4.5
STF-2/16	16	M4.5
STF-3/06	6.0	M6
STF-3/08	8.0	M6
STF-3/10	10	M6
STF-3/11.5	11.5	M6
STF-3/13	13	M6
STF-3/16	16	M6



DRILL EXTENSION

Equipment made by stainless steel manufactured with two connection types.

Article	Diameter	Total Length	Internal Cooling	Connection
AIP-1	2.35/4.5	30	optional	Bayonet
AIP-3	2.35/4	30	optional	Easy



HAND TISSUE PUNCH

Tissue Punch made by stainless steel apt to create a circular incision of the gingival mucosa.

Available with or without probe attachment.

Article	Diameter	Total Length	Probe
MUC-300	Ø3.0	32	No
MUC-350	Ø3.5	32	No
MUC-410	Ø4.1	32	No
MUC-450	Ø4.5	32	No
MUC-510	Ø5.1	32	No
MUC-300C	Ø3.0	32	Yes
MUC-350C	Ø3.5	32	Yes
MUC-410C	Ø4.1	32	Yes
MUC-450C	Ø4.5	32	Yes
MUC-510C	Ø5.1	32	Yes



BONE EXPANDER

Device made by stainless steel which allows to expand and compact bone. Consisting of 4 pieces of different diameters.

Article	External Connection / Internal Connection
KIT-OST1	4mm Square / 2mm Hexagon
KIT-OST2	4mm Square / 2.5mm Hexagon
KIT-OST3	4mm Square with O-ring



CONTRA-ANGLE SCREW DRIVER

Screw driver made by stainless steel available in different dimensions.

Article	Connection	Total Length
ACE-1	1.20mm Hexagon	22
ACE-2	1.25mm Hexagon	22
ACE-3	2.00mm Hexagon	22
ACE-4	0.90mm Hexagon	22
ACE-5	2.50mm Hexagon	22
ACE-7	"Special Star"	22
ACE-8	Internal 4.00mm Square	23
ACE-9	2.40mm Hexagon with o-ring	24
ACE-10	1.25mm Square	22
ACE-11	1.25mm Square	27
ACE-12	2.40mm Hexagon	22
ACE-1L	1.20mm Hexagon	27
ACE-2L	1.25mm Hexagon	27
ACE-4L	0.90mm Hexagon	27
ACE-12L	2.40mm Hexagon	27



HAND SCREW DRIVER

Hand screw driver made by stainless steel available in different dimensions.

Article	Connection	Total Length
AME-1	1.25mm Hexagon	20
AME-2	1.25 mm Hexagon	25
AME-5	1.20 mm Hexagon	20
AME-6	1.20 mm Hexagon	25
AME-7	1.50 mm Hexagon	20
AME-8	1.50 mm Hexagon	25
AME-9	1.80 mm Hexagon	20
AME-10	1.80 mm Hexagon	25
AME-11	2.50 mm Hexagon	20
AME-12	2.50 mm Hexagon	25
AME-13	0.90 mm Hexagon	20
AME-14	0.90 mm Hexagon	25
AMQ-1	1.00mm Square	20
AMQ-2	1.00mm Square	25



CONTRA-ANGLE SCREW DRIVER HOLDER

Equipment made by stainless steel which allows to manually use every devices with handpiece connection.

Article	Diameter	Total Length
PCM-1	Ø10	13



TORQUE WRENCH ADAPTER

Equipment made by stainless steel which allows to use every devices with handpiece connection by a torque wrench and ratchet.

Article	Torque Wrench Connection	Device Connection
ADT-1	Ø7 with small wheel	Handpieces
ADT-3	Ø7mm open	Handpieces
ADT-4	Ø7	4mm Square
ADT-5	4mm Square open	Handpieces
ADT-6	6mm Hexagon open	Handpieces
ADT-7	6.35mm Hexagon open	Handpieces
ADT-8	6.35mm Hexagon	4mm Square

We manufactured all types of connection on demand.



RATCHER WRENCH

Device made by stainless steel available in different connection types.

Article	Connection
CRF-1	4.00mm Square
CRF-2	Ø7.00mm
CRF-3	3.00mm Square
CRF-4	3.55mm Hexagon
CRF-5	6.00 mm Hexagon
CRF-6	3.80 mm Hexagon
CRF-7	4.65 mm Hexagon
CRF-8	3.00 mm Hexagon
CRF-11	6.35 mm Hexagon



TORQUE WRENCH

Device made by stainless steel available in different types of torques and connections.

Article	Connection	Torque
CDM-1	4.00mm Square	10-35Ncm
CDM-2	Ø7.00mm	10-35Ncm
CDM-3	3.00mm Square	10-35Ncm
CDM-4	3.55mm Hexagon	10-35Ncm
CDM-5	6.00mm Hexagon	10-35Ncm
CDM-6	4.00mm Square	25-55Ncm
CDM-7	Ø7.00mm	25-55Ncm
CDM-8	3.00mm Square	25-55Ncm
CDM-9	3.55mm Hexagon	25-55Ncm
CDM-10	6.00mm Hexagon	25-55Ncm
CDM-11	6.35mm Hexagon	10-35Ncm
CDM-12	6.35mm Hexagon	25-55Ncm
CDM-13	3.55mm Hexagon	25-55Ncm
CDM-14	4.00mm Hexagon	25-55Ncm
CDM-15	4.00mm Square	40-70Ncm
CDM-16	3.00mm Hexagon	10-35Ncm
CDM-17	3.00mm Hexagon	40-70Ncm
CDM-18	4.00mm Square	30-60Ncm
CDM-19	3.70mm Hexagon	25-55Ncm
CDM-20	4.00mm Square	5-10Ncm
CDM-21	4.70mm Hexagon	25-55Ncm
CDM-22	6.35mm Hexagon	20-45Ncm
CDM-23	6.00mm Hexagon	25-60Ncm
CDM-24	6.35mm Hexagon	10-45Ncm
CDM-25	4.00mm Square	10-45Ncm
CDM-26	3.80mm Hexagon	25-75Ncm
CDM-27	Ø7.00mm	10-40Ncm
CDM-101	4.00mm Square	40-110Ncm

We manufactured all types of connection on demand.



SURGICAL BOX

There are several types of surgical boxes for every needs, with the possibility of customization.

